



The right probe for any application

The supreme 3D/4D ultrasound, ACCUVIX V20 optional probes provide exceptional imaging performance in all applications. Advanced technology of automatic 3D volume reconstructing and expanded functional bandwidth delivers increased amount of clinical information and qualified diagnosis.



www.medison.com | marketing@medison.com



see it all
MEDISON

PG-MD-V20-ICM-06.30.2008



















Probe Guide

Version 1.0

see it all
MEDISON

Probe Guide

																
TRANSDUCER	C2-5EL	C2-6IC	C3-7IM	EC4-9IS	L5-13IS	L7-16IS	L5-12/50EP	L4-7EL	3D2-6ET	3D4-8ET	3D4-9ES	3D5-9EK	P2-4AC	P3-5AC	2.0CW	4.0CW
Array Type	Curved	Curved	Curved	Curved	Linear	Linear	Linear	Linear	Curved	Curved	Curved	Curved	Phased	Phased	CW	CW
Sector Angle/ Field of View	76°	58°	78°	149°	38mm	38mm	50mm	38mm	84°	84°	128°	128°	19mm	14mm	N/A	N/A
Radius of Curvature	40mm	60mm	50mm	10mm	Flat	Flat	Flat	Flat	40mm	40mm	10mm	10mm	Flat	Flat	N/A	N/A
Broadband Frequency Range	2MHz~5MHz	2MHz~6MHz	3MHz~7MHz	4MHz~9MHz	5MHz~13MHz	7MHz~16MHz	5MHz~12MHz	4MHz~7MHz	2MHz~6MHz	4MHz~8MHz	4MHz~9MHz	5MHz~9MHz	2MHz~4MHz	3MHz~5MHz	2MHz	4MHz
APPLICATION																
Abdomen	●	●	●						●	●						
Gynecology	●	●	●	●					●	●	●	●				
OB	●	●	●	●					●	●	●	●				
OB Early	●	●	●	●					●	●						
Fetal Heart	●	●	●	●					●	●						
Musculoskeletal					●	●	●	●								
Urology	●	●	●	●							●	●				
Pediatric					●	●	●	●					●	●		
Small Part					●	●	●	●								
Breast					●	●	●	●								
Vascular					●	●	●	●								
Cardiology													●	●	●	●
TCD													●	●		
Neonatal													●	●		
NEEDLE GUIDE	DISPOSABLE	DISPOSABLE	REUSABLE	REUSABLE	DISPOSABLE	DISPOSABLE	DISPOSABLE	DISPOSABLE	DISPOSABLE	DISPOSABLE	REUSABLE	REUSABLE	N/A	N/A	N/A	N/A
Biopsy Kit	BP-KIT-010	BP-KIT-009	BP-C3-7IM	BP-KIT-003	BP-KIT-012	BP-KIT-012	BP-KIT-020	BP-KIT-014	BP-KIT-017	BP-KIT-017	VF-BP-PEC47	VF-BP-PEC47				
	